# Levenhuk MED 1000 Series Biological Microscopes

Levenhuk MED 1000B / 1000T / D1000T

- EN User Manual
- BG Ръководство за потребителя
- CZ Návod k použití
- DE Bedienungsanleitung
- ES Guía del usuario
- HU Használati útmutató
- IT Guida all'utilizzo
- PL Instrukcja obsługi
- RU Инструкция по эксплуатации

*Zoom&Joy*

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## General Information

Levenhuk MED 1000 Series biological microscopes are safe for health, life and property of the consumer and the environment when properly used, and meets the requirements of international standards. These microscopes are designed for observing transparent objects in the transmitted light using the bright field method. They perform well in clinical research and tests; teaching demonstrations; bacterioscopy and cytoscopy in medical and health establishments, laboratories, universities, and schools; and may be used for scientific research in agriculture and microbiology.

## The kit includes

- microscope body
- eyepiece head: binocular (Levenhuk MED 1000B), trinocular (Levenhuk MED 1000T and D1000T)
- infinity-corrected planar objectives: 4x, 10x, 40xs, 100xs (oil)
- eyepieces: WF10x/22mm (2pcs)
- blue filter
- vial of immersion oil
- C-mount for camera (Levenhuk MED 1000T and D1000T)
- dust cover
- power cable
- LED lamp
- user manual and warranty

The Levenhuk MED D1000T kit additionally includes: 14Mpx digital camera, a USB cable and a CD with software & drivers.

## Optional accessories

You can increase the possibilities of your microscope by equipping it with the optional accessories of the Levenhuk MED series: eyepieces of different magnification, reticle eyepieces, planar objectives, dark-field condenser, LCD microscope camera, etc. Visit our website to choose and buy optional accessories: [www.levenhuk.com](http://www.levenhuk.com)

## Microscope assembly

1. **Eyepiece head.** Loosen the locking screw on the bottom of the microscope body. Install the eyepiece head and tighten the locking screw.
2. **Eyepiece.** There is a special ring to fix the eyepiece firmly in the tube of eyepiece head. Make sure to install the proper end of the eyepiece in the tube.
3. **Installing objective lenses.** To install the objective lenses, screw them in by turning clockwise. Attach objective lenses to the revolving nosepiece in ascending order of magnification.
4. **Condenser.** The condenser was mounted, focused and centered at the factory. There is no need to make any adjustments.
5. **Digital camera** (Levenhuk MED D1000T). Insert the camera into the C-mount and tighten the locking screw. Then, insert the C-mount into the trinocular tube.

## Use

- **Illumination.** Connect your microscope to a power outlet. Turn the power on; the power indicator will come on. You can change the level of brightness with the brightness adjustment wheel.
- **Iris diaphragm** is incorporated into the condenser. Adjust the iris diaphragm to match the numerical aperture of different objectives. Open or close the iris diaphragm with the adjusting knob. When the knob is on the right, the iris diaphragm is fully opened.
- **Placing the specimen.** Place the specimen on the stage and cover it with glass. Make sure to fix the specimen securely with the stage clips. Move the stage in horizontal and vertical directions to point the objective lens at the specimen.
- **Focusing.** Select the objective lens with the lowest magnification rotating the revolving nosepiece. Raise the stage turning the coarse focus knob. While observing through the objective lens, turn the coarse focus knob slowly to let the stage down until you see an image of the observed specimen. Turning the fine focus knob, make the image sharp. **CAUTION:** The objective should not touch the specimen; otherwise, the objective and/or the specimen might be damaged.
- **Adjusting the eyepiece head.** Adjust the eyepiece tubes to make the view in the left and right eyepiece even. Start the diopter adjustment with “0” value and continue until the image becomes sharp.
- **Immersion oil.** For observing at high magnifications with a 100x objective lens, you have to use immersion oil. Place a drop of immersion oil on the specimen. Turn the revolving nosepiece to place the 100x objective lens in its working position. Raise the microscope stage with a coarse focus knob until the bottom of the objective lens touches the oil drop. Move the revolving nosepiece back and forth to get rid of air bubbles in the oil. Then, fix the objective lens in its working position. Make sure that oil fills the entire space between the objective lens and the specimen.
- **Changing the bulb.** The microscope is equipped with built-in LED illumination. Change the bulb only if necessary.

## Digital camera

Levenhuk MED D1000T comes with a 14Mpx digital camera. The camera allows you to observe specimens in fine detail and true colors on your PC monitor and save images on the hard drive. Connect the camera to your PC via a USB cable that also acts as a power cable.

## Digital camera specifications

Max. resolution (still images)	4096x3288px
Megapixels	14
Sensor	1/2.3" CMOS
Mounting location	third eyepiece (with a C-mount)
Pixel size	1.4x1.4µm
Sensitivity	0.724V/lux.sec@550nm
Spectral range	380–650nm
Image format	*.jpg, *.bmp, *.png, *.tif
Video format	recording: *.wmv, *.avi, *h.264 (Win 8 and above), *h.265 (Win 10 and above)
Exposure	ERS (Electronic Rolling Shutter)
White balance	auto/manual
Exposure control	auto/manual
Software features	image size, brightness, exposure control
Port	USB 2.0, 480Mb/s
System requirements	Mac OS 10.12, Linux Ubuntu 14.04, Windows XP/Vista/7/8/10 (32 and 64 bits), CPU 2.8GHz Intel Core 2 or higher, RAM not less than 2Gb, USB port 2.0, CD-Rom
Software	Levenhuk
Camera power supply	USB 2.0 cable

## Microscopes specifications

	Levenhuk MED 1000B	Levenhuk MED 1000T	Levenhuk MED 1000T
Type	biological	biological	digital, biological
Eyepiece head	binocular, 360° rotatable, inclined at 30°	trinocular, 360° rotatable, inclined at 30°	
Revolving nosepiece		quadruple	
Magnification		40–1000x	
Objectives	infinity-corrected planar objectives 4x, 10x, 40xs, 100xs (oil)		
Eyepieces	WF10x/22mm (2pcs)		
Optics material	glass		
Body material	aluminum, plastic		
Eyepiece tube diameter	30mm	30mm (binocular eyepiece head) 23.2mm (third vertical tube)	
Stage	double layer mechanical stage, 216x150mm, with slide holder and graphite-coated surface		
Stage movement range	75mm/55mm		
Condenser	Abbe N.A. 1.25 iris diaphragm and filter		
Focus system	coaxial, coarse (37.7mm) and fine (0.001mm)		
Illumination	Koehler illumination lower (3W LED, brightness adjustment)		
Power source	AC adapter 100–220V / 50–60Hz		
Camera	–	–	14Mpx

Levenhuk reserves the right to modify or discontinue any product without prior notice.

## Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

## Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase.

Levenhuk will repair or replace such product or part thereof which, upon inspection by Levenhuk, is found to be defective in materials or workmanship. As a condition to the obligation of Levenhuk to repair or replace such product, the product must be returned to Levenhuk together with proof of purchase satisfactory to Levenhuk.

This warranty does not cover consumable parts, such as bulbs (electrical, LED, halogen, energy-saving and other types of lamps), batteries (rechargeable and non-rechargeable), electrical consumables etc. For further details, please visit our web site: [www.levenhuk.com/warranty](http://www.levenhuk.com/warranty)

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2018 Levenhuk, Inc. All rights reserved.
Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA